



Software Release Notes for nordicICE v2.3.1

Date: September 11th, 2008
Approved by: Anders Strand Vestbø

NOTE: The most current documentation for released products is available on <http://www.nordicimaginglab.com>. The online documentation may contain updated information on modifications made and issues reported after this document was issued.

Introduction

These release notes describe the bug-fixes and enhancements included in nordicICE, software version v2.3.1 released on September 11th, 2008.

Overview

This update basically implements numerous minor bug-fixes within the DCE/Permeability module.

It also provides an initial beta/test-version of the Report Module. This module can be installed by choosing the 'Custom' installation scheme. Note that if you are updating nordicICE from v2.3.0, you need to run the installer package a second time after updating in order to be able to install the module. Note also that no Help manual is included yet for this module.

Supplementary Notes

The software can be downloaded and installed from the company website:

<http://www.nordicimaginglab.com>

Contact

For questions or problems please contact support:

support@nordicimaginglab.com